




# DynaMesh<sup>®</sup>

Implants for the surgical treatment of  
Female Urinary Incontinence  
Female Pelvic Organ Prolapse  
Male Urinary Incontinence



■ made  
■ in  
■ Germany

Tailored Implants  
Made of **PVDF**



## Profile and Competences



FEG Textiltechnik Forschungs- und Entwicklungsgesellschaft mbH (FEG) was established in Aachen in 1992 and since then has gained a reputation for exceptional innovations in the field of textile technology.

The company's dedicated and highly qualified staff, in-house research and development capacities as well as intelligent production facilities combine to make FEG Textiltechnik today's leading manufacturer of textile surgical implants in Germany. Under the brand name **DynaMesh®**, an internationally protected trademark, FEG's award-winning implants are successfully marketed in numerous countries around the world.

Constant and close contact with major scientific, medical and technical institutions ensures that FEG's high-quality products meet the latest requirements in terms of patient comfort and surgical handling. The sophisticated quality management system at FEG Textiltechnik is fully certified to DIN EN ISO 13485 for the manufacture of medical devices. All of FEG's products are CE approved (CE 0123) and are approved under relevant national regulations.

Focusing on its core competences and expertise in textile implants, FEG Textiltechnik will continue to set technical benchmarks in the future.

# DynaMesh®

# Milestones

1992 Founding of:



2003 Certification of:



2011 Development of MRI-visible technology



2014 New 4,200m<sup>2</sup> offices & production plant



2020 Additional 600m<sup>2</sup> production/storage capacity



1994 Active in medical technology

2004 First implant for the surgical treatment of hernias

2005 First implant for the surgical treatment of female urinary incontinence

2006 First implant for the surgical treatment of female pelvic organ prolapse

2007 First implant for the surgical treatment of parastomal hernias

2008 First implant for the surgical treatment of male urinary incontinence




2016 First implant for the surgical treatment of hiatal hernias

### Business Fields:

- Implants for the surgical treatment of:
  - Hernias
  - Female Pelvic Organ Prolapse
  - Female Urinary Incontinence
  - Male Urinary Incontinence
- Sales in over 50 countries
- More than 70 employees

# Implants 'made in Germany'

## DynaMesh<sup>®</sup>



**Filament**

**Warp-Knitted Fabric**

**Implant**

Spinning

Warp-Knitting


Finishing

Washing, Sterilisation

From Thread to Implant:

Full quality control along the entire production chain.

Development and manufacturing in Aachen, Germany




# Biocompatibility

PVDF mesh structures have good biocompatibility (assessed according to ISO 10993) and show significantly lower granuloma formation (scar tissue) [1<sup>A</sup>,2<sup>A</sup>,4<sup>A</sup>,68<sup>A</sup>,100<sup>B</sup>]. Therefore, the risk of undesirable foreign body reactions is minimised.

## Cross-Sectional View

A comparison of different granuloma thicknesses





[#] Reference "#" (see "References")  
 [#<sup>A</sup>] Reference "#" (see "References"),  
 "A": limitation "animal trial"  
 [#<sup>B</sup>] Reference "#" (see "References"),  
 "B": limitation "in-vitro trial"

Vi001 de	DynaMesh® Implantate - Animation: Fremdkörperreaktion - Vergleich zwischen PVDF und PP <a href="https://de.dyna-mesh.com/Vi001de">https://de.dyna-mesh.com/Vi001de</a>	
Vi001 en	DynaMesh® Implants - Animation: Foreign Body Reaction - Comparison of PVDF and PP <a href="https://de.dyna-mesh.com/Vi001en">https://de.dyna-mesh.com/Vi001en</a>	

# Ageing Resistance

PVDF has been used as a surgical suture material for many decades with great success, even in the most demanding areas of application such as ophthalmology and cardiology [91].

Long-term data with observation periods of up to seven years prove that: The condition of the PVDF surface remains unchanged, filaments are still stable, nothing becomes brittle [101,2<sup>A</sup>,5<sup>B</sup>,27<sup>A</sup>,52<sup>B</sup>,93<sup>A</sup>].



**PVDF-polyvinylidene fluoride**





**PP-polypropylene**

[#] Reference "#" (see "References")  
 [#<sup>A</sup>] Reference "#" (see "References"), "A": limitation "animal trial"  
 [#<sup>B</sup>] Reference "#" (see "References"), "B": limitation "in-vitro trial"

## Dynamometry

Textile implants must reinforce tissue. They have to cushion different forces - including the extreme stresses associated with coughing, sneezing and laughing. What is needed therefore, is a good interaction between stability and elasticity.


### The behaviour of abdominal walls with different mesh implants under load



A study of explanted abdominal walls (source: Aachen University Hospital, Germany)

## Tear Propagation Resistance


The multiple meshing technique in warp-knitted\* DynaMesh® structures minimises the risk of the zipper effect (once torn, the structure tears further).  
The load-specific adapted tear resistance is one of the key properties of DynaMesh® implants.



DynaMesh® products are not woven or conventionally knitted, but warp-knitted\*. This technology, unlike any other, makes it possible to make specific variations in the shape and structure of a textile implant, which means that we can construct features with different characteristics in different places within the structure. It is impossible to achieve a more accurate adaptation of implants to the relevant indication.

### \*Warp-Knitted Fabric


Warp-knitted fabrics are a type of knitted fabric. They are produced industrially on warp-knitting machines via stitch formation from thread systems.






## Specially Warp-knitted Selvedges


DynaMesh® products are not simply cut from a flat piece of mesh.  
Special warp-knitting machines produce smooth selvedges  
(no 'sawtooth' edges!).



Conventional Mesh




**DynaMesh®-LICHTENSTEIN**



**DynaMesh®-SIS soft**

Also in three-dimensionally  
shaped implants




**DynaMesh®-IPST-R visible**

## Effective Porosity

During incorporation, the filaments are enclosed by an internal and external granuloma. When filament distance is too small, there is a risk that the whole intervening space will be filled with scar tissue (closed pores). Sufficiently large pores can prevent this [8C].

**How is this prevented?** PP implants must have a **pore diameter of at least 1 mm** in all directions, even under load! In the case of **PVDF** implants, a diameter of **0.6 mm** is already sufficient to keep the pores open due to the lower granuloma thickness. Only in this way can local autochthonous tissue form through a pore [6,68<sup>A</sup>,105<sup>A</sup>].



[#] Reference "#" (see "References")  
[<sup>A</sup>] Reference "#" (see "References"), "<sup>A</sup>": limitation "animal trial"  
[<sup>C</sup>] Reference "#" (see "References"), "<sup>C</sup>": limitation "results based on the analysis of explants"

**Textile porosity** refers to the permeable component of a mesh implant **before** the body has reacted to it.

**Effective porosity** refers to the permeable component of a mesh implant **after** the body has reacted to it.

**Rule of thumb:**


A 'pore' less than 1 mm in diameter is closed by the body with scar tissue  
= 0% effective porosity

## DynaMesh® visible

Conventional mesh implants are mostly undetectable in diagnostic radiology.

DynaMesh® visible implants can be detected using magnetic resonance imaging (MRI) [7,29<sup>A</sup>,51,54,56,62,69-71,76,90] - both in standard sequences and in high-resolution, three-dimensional images and even films. As such, position and condition of the implant can be determined reliably and accurately.

DynaMesh® visible is the world's first technology to visualise textile implants. The PVDF filament is mixed with ferromagnetic micro-pigments according to a proprietary process that ensures optimum pigment incorporation. This innovation has won an award from the German Federal Ministry of Education and Research (BMBF 01EZ 0849).



Award-winner in the innovation competition hosted by the




The development was sponsored by the  
German Federal Ministry  
of Education and Research  
(BMBF 01EZ 0849)

[#] Reference "#" (see "References")


[#<sup>A</sup>] Reference "#" (see "References"), "A": limitation "animal trial"

# DynaMesh® visible




## DynaMesh®-PRS visible 3-dimensional remodelling





## DynaMesh®-ENDOLAP visible 3-dimensional remodelling



During MRI scans, the part of the body being analysed is scanned step-by-step and deconstructed into many 'wafer-thin optical slices'. At the end, these 'slices' are reconstructed to form 3-dimensional images or motion sequences (remodelling).


Vi069XX	DynaMesh® MRI - Animation: MRI Reconstruction with DynaMesh®-PRP visible <a href="https://de.dyna-mesh.com/Vi069xx">https://de.dyna-mesh.com/Vi069xx</a>	
Vi067XX	DynaMesh® MRI - Animation: MRI Reconstruction with DynaMesh®-PRS visible <a href="https://de.dyna-mesh.com/Vi067xx">https://de.dyna-mesh.com/Vi067xx</a>	
Vi032XX	DynaMesh®-ENDOLAP visible - Animation: MRI visible - 3D Implant Remodelling <a href="https://de.dyna-mesh.com/Vi032xx">https://de.dyna-mesh.com/Vi032xx</a>	

**LiSTO.academy** is a comprehensive platform providing excellent customised surgical education and trainings. Through cooperations with experienced surgeons and renowned experts from around the world, **LiSTO.academy** enables physicians to achieve best patient outcomes using DynaMesh® implants.



Visit us online:  
[www.listo.academy](http://www.listo.academy)






DynaMesh®-SIS implants are designed as a midurethral sling for soft tissue reinforcement of the pelvic floor as part of the surgical treatment of stress urinary incontinence caused by a hypermobile urethra and/or intrinsic sphincter deficiency.

## DynaMesh®-SIS


<b>DynaMesh®-SIS</b>	01 cm x 50 cm	PV211056F1	BX = 1 piece
		PV211056F3	BX = 3 pieces
<b>DynaMesh®-SIS soft</b>	01 cm x 50 cm	PV411056F1	BX = 1 piece
		PV411056F3	BX = 3 pieces
<b>DynaMesh®-SIS visible</b>	01 cm x 50 cm	PV471056F1	BX = 1 piece
		PV471056F3	BX = 3 pieces




### Use and Properties

Product	DynaMesh®-SIS	DynaMesh®-SIS soft	DynaMesh®-SIS visible
Field of application	stress urinary incontinence (SUI)		
Surgical access	transvaginal		
Surgical technique	TVT - retropubic - inside-out / TOT - transobturator - inside-out - outside-in		
Fixation	none		
Specially Warp-knitted Selvedges		●	
Shape stability [TR1,TR12]		●	
Defined elasticity [TR10]		●	
Visible technology	●	●	●
Polymer (monofilament)		PVDF	
Biocompatibility [1 <sup>A</sup> ,2 <sup>A</sup> ,4 <sup>A</sup> ,68 <sup>A</sup> ,100 <sup>A</sup> ,TR1]		●	
Ageing resistance [101,2 <sup>A</sup> ,5 <sup>B</sup> ,52 <sup>B</sup> ,93 <sup>A</sup> ,27 <sup>A</sup> ]		●	
Classification (Klinge's classification [8]) [TR1 1]		1a	

● Applies to all product sizes  
● Does not apply  
[#] Reference "#" (see "References")  
[#<sup>A</sup>] Reference "#" (see "References"), "A": limitation "animal trial"  
[#<sup>B</sup>] Reference "#" (see "References"), "B": limitation "in-vitro trial"  
[TR#] Internal test-report (see "internal test-report references")



**Retropubic**  
(inside-out / bottom-up)



**Transobturator**  
(inside-out & outside-in)



**DynaMesh®-SIS** implants are positioned using the inside-out technique in case of a retropubic tape position, and using the outside-in or inside-out technique in case of a transobturator tape position.

**DynaMesh®-SIS** implants have a thread on both ends of the sling, which aids fixation to the surgical instrument.

Several reusable instruments are available separately to assist the positioning of **DynaMesh®-SIS** implants:



**DynaMesh®-ISR01:**

Instrument for retropubic positioning of **DynaMesh®-SIS** implants through transvaginal access using the inside-out technique.



**DynaMesh®-IST01/-IST02/-IST03:**

Instrument set consisting of two instruments (right and left side) for transobturator positioning through transvaginal access using the inside-out or outside-in technique.

Diameter: 5 - 7 cm



**DynaMesh®-IVT01:**


Instrument for transobturator positioning through transvaginal access using the outside-in technique.

VI040xx DynaMesh®-SIS - Animation: Surgical Treatment of Stress Urinary Incontinence - SUI - TVT 8/4  
<https://de.dyna-mesh.com/Vi040xx>



Distributed by:


**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: info@zisanamed.ch



DynaMesh®-SIS direct implants are designed as a midurethral sling for soft tissue reinforcement of the pelvic floor as part of the surgical treatment of stress urinary incontinence caused by a hypermobile urethra and/or intrinsic sphincter deficiency.

## DynaMesh®-SIS direct

DynaMesh®-SIS direct	01 cm x 50 cm	PV211050F1	BX = 1 piece
		PV211050F3	BX = 3 pieces
DynaMesh®-SIS direct soft	01 cm x 50 cm	PV411050F1	BX = 1 piece
		PV411050F3	BX = 3 pieces
DynaMesh®-SIS direct visible	01 cm x 50 cm	PV471050F1	BX = 1 piece
		PV471050F3	BX = 3 pieces




### Use and Properties

Product	DynaMesh®-SIS direct	DynaMesh®-SIS direct soft	DynaMesh®-SIS direct visible
Field of application	stress urinary incontinence (SUI)		
Surgical access	transvaginal		
Surgical technique	TOT - transobturator - outside-in		
Fixation	none		
Specially Warp-knitted Selvedges		●	
Shape stability [TR1,TR12]		●	
Defined elasticity [TR10]		●	
Visible technology	●	●	●
Polymer (monofilament)		PVDF	
Biocompatibility [1 <sup>A</sup> ,2 <sup>A</sup> ,4 <sup>A</sup> ,68 <sup>A</sup> ,100 <sup>A</sup> ,TR1]		●	
Ageing resistance [101,2 <sup>A</sup> ,5 <sup>B</sup> ,52 <sup>B</sup> ,93 <sup>A</sup> ,27 <sup>A</sup> ]		●	
Classification (Klinge's classification [8]) [TR11]		1a	


● Applies to all product sizes  
● Does not apply  
[#] Reference "#" (see "References")  
[#<sup>A</sup>] Reference "#" (see "References"), "A": limitation "animal trial"  
[#<sup>B</sup>] Reference "#" (see "References"), "B": limitation "in-vitro trial"  
[TR#] Internal test-report (see "internal test-report references")






DynaMesh®-SIS direct   DynaMesh®-SIS direct soft   DynaMesh®-SIS direct visible

Transobturator (outside-in)



DynaMesh®-SIS direct implants are positioned using the outside-in technique in a transobturator tape position.

Several reusable instruments are available separately to assist the positioning of DynaMesh®-SIS direct implants:



Diameter: 5 - 7 cm

**DynaMesh®-IST01/-IST02/-IST03:**

Instrument set consisting of two instruments (right and left side) for transobturator positioning through transvaginal access using the outside-in technique.



**DynaMesh®-IVT01:**


Instrument for transobturator positioning through transvaginal access using the outside-in technique.

VI045en	DynaMesh®-SIS direct - Animation: SUI Treatment - Transobturator (out/in) - TOT <a href="https://de.dyna-mesh.com/Vi045en">https://de.dyna-mesh.com/Vi045en</a>	
VI021xx	DynaMesh®-SIS direct - Animation: SUI Treatment - Transobturator (out/in) - TOT 8/4 <a href="https://de.dyna-mesh.com/Vi021xx">https://de.dyna-mesh.com/Vi021xx</a>	

Distributed by:

**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: info@zisanamed.ch


Female Pelvic Organ Prolapse  
Vaginal/Cervical Stump Prolapse  
or Uterine Prolapse



DynaMesh®-PR implants have been specially developed for pelvic floor reconstruction, in laparoscopic or open surgical technique, and serve to support and stabilise fascial structures, connective tissue and ligaments.

## DynaMesh®-PR

<b>DynaMesh®-PR soft</b>	04 cm x 23 cm	PV500423F1	BX = 1 piece
		PV500423F3	BX = 3 pieces
<b>DynaMesh®-PR visible</b>	04 cm x 23 cm	PV700423F1	BX = 1 piece




### Use and Properties

Product	DynaMesh®-PR soft	DynaMesh®-PR visible
Field of application	vaginal/cervical stump or uterine prolapse, concomitant cystocele/rectocele	
Surgical access	laparoscopic / open	
Surgical technique	colposacropexy / cervicosacropexy / hysterovacropexy unilateral	
Fixation on vagina / cervix	sutures	
Fixation on sacrum	sutures / tacks	
Specially Warp-knitted Selvedges		●
Shape stability		●
Defined elasticity		●
Visible technology	●	●
Polymer (monofilament)		PVDF
Biocompatibility		●
Ageing resistance		●
Dynamometry		●
Tear propagation resistance		●
Classification (Klinge's classification [8])		1 a

● Applies to all product sizes  
● Does not apply


DynaMesh®-PR implants are used in the surgical treatment of the vaginal/cervical stump or uterine prolapse, as well as in the treatment of a concomitant cystocele/rectocele.

**Application Examples:**




Colpo-/cervicosacropexy

- unilateral
- fixation on vaginal/cervical stump




Colpo-/cervicosacropexy

- unilateral
- fixation on vaginal/cervical stump and anterior mesh plasty for concomitant cystocele




Colpo-/cervicosacropexy

- unilateral
- fixation on vaginal/cervical stump and posterior mesh plasty for concomitant rectocele




Hysterosacropexy

- unilateral
- posterior cervical fixation



Hysterosacropexy


- unilateral
- posterior cervical fixation and posterior mesh plasty for concomitant rectocele

VI086xx	DynaMesh®-PR - Animation: Colposacropexy <a href="https://de.dyna-mesh.com/Vi086xx">https://de.dyna-mesh.com/Vi086xx</a>	
---------	--	---

Distributed by:

**Zi+Sanamed GmbH**  
 Bernstrasse 18  
 CH-2555 Brügg  
 phone +41 793310767  
 email: info@zisanamed.ch


Female Pelvic Organ Prolapse  
Vaginal/Cervical Stump Prolapse  
or Uterine Prolapse



DynaMesh®-PRR implants have been specially developed for pelvic floor reconstruction, in laparoscopic or open surgical technique, and serve to support and stabilise fascial structures, connective tissue and ligaments.

## DynaMesh®-PRR

DynaMesh®-PRR soft	02/04 cm x 23 cm	PV360423F1	BX = 1 piece
		PV360423F3	BX = 3 pieces
DynaMesh®-PRR visible	02/04 cm x 23 cm	PV760423F1	BX = 1 piece




### Use and Properties

Product	DynaMesh®-PRR soft	DynaMesh®-PRR visible
Field of application	vaginal/cervical stump or uterine prolapse, concomitant cystocele/rectocele	
Surgical access	laparoscopic / open	
Surgical technique	colposacropexy / cervicosacropexy / hysterocropexy unilateral	
Fixation on vagina / cervix	sutures	
Fixation on sacrum	sutures / tacks	
Specially Warp-knitted Selvedges		●
Shape stability		●
Defined elasticity		●
Visible technology	●	●
Polymer (monofilament)		PVDF
Biocompatibility		●
Ageing resistance		●
Dynamometry		●
Tear propagation resistance		●
Classification (Klinge's classification [8])		1a

● Applies to all product sizes  
● Does not apply


DynaMesh®-PRR implants are used in the surgical treatment of the vaginal/cervical stump or uterine prolapse, as well as in the treatment of a concomitant cystocele/rectocele.

**Application Examples:**




Colpo-/cervicosacropexy

- unilateral
- fixation on vaginal/cervical stump




Colpo-/cervicosacropexy

- unilateral
- fixation on vaginal/cervical stump and anterior mesh plasty for concomitant cystocele




Colpo-/cervicosacropexy

- unilateral
- fixation on vaginal/cervical stump and posterior mesh plasty for concomitant rectocele



Hysteriosacropexy

- unilateral
- posterior cervical fixation




Hysteriosacropexy

- unilateral
- posterior cervical fixation and posterior mesh plasty for concomitant rectocele

VI083xx	DynaMesh®-PRR - Animation: Colposacropexy <a href="https://de.dyna-mesh.com/Vi083xx">https://de.dyna-mesh.com/Vi083xx</a>	
Vi062xx	DynaMesh®-PRR - Animation: Hysteriosacropexy <a href="https://de.dyna-mesh.com/Vi062xx">https://de.dyna-mesh.com/Vi062xx</a>	






Distributed by:  
  
**Zi+Sanamed GmbH**  
 Bernstrasse 18  
 CH-2555 Brügg  
 phone +41 793310767  
 email: info@zisanamed.ch



**DynaMesh®-PRS** implants have been specially developed for pelvic floor reconstruction, in laparoscopic or open surgical technique, and serve to support and stabilise fascial structures, connective tissue and ligaments. The implants are used in the surgical treatment of the vaginal or cervical stump prolapse, as well as in the treatment of a concomitant cystocele and/or rectocele.

## DynaMesh®-PRS

<b>DynaMesh®-PRS soft</b>	<sup>1)</sup> 02 cm x 16 cm	PV350216F1	BX = 1 piece
<b>DynaMesh®-PRS soft</b>	<sup>2)</sup> 03 cm x 23 cm	PV350323F1	BX = 1 piece
<b>DynaMesh®-PRS visible</b>	<sup>3)</sup> 03 cm x 23 cm	PV750323F1	BX = 1 piece
<b>DynaMesh®-PRS visible</b>	<sup>4)</sup> 3.3 cm x 24 cm	PV750424F1	BX = 1 piece
		PV750424F10	BX = 10 pieces
<b>DynaMesh®-PRS visible</b>	<sup>5)</sup> 04 cm x 20 cm	PV750420F1	BX = 1 piece
		PV750420F10	BX = 10 pieces
<b>DynaMesh®-PRS soft</b>	<sup>6)</sup> 05 cm x 27 cm	PV350527F1	BX = 1 piece




VI046xx	DynaMesh®-PRS - Animation: Colposacropexy <a href="https://de.dyna-mesh.com/Vi046xx">https://de.dyna-mesh.com/Vi046xx</a>	
VI048xx	DynaMesh®-PRS - Animation: Colposacropexy <a href="https://de.dyna-mesh.com/Vi048xx">https://de.dyna-mesh.com/Vi048xx</a>	
VI067xx	DynaMesh® MRI - Animation: MRI Reconstruction with DynaMesh®-PRS visible <a href="https://de.dyna-mesh.com/Vi067xx">https://de.dyna-mesh.com/Vi067xx</a>	

## Use and Properties


Product	DynaMesh®-PRS soft 02 cm x 16 cm <sup>(1)</sup> 03 cm x 23 cm <sup>(2)</sup> 05 cm x 27 cm <sup>(6)</sup>	DynaMesh®-PRS visible 03 cm x 23 cm <sup>(3)</sup> 3.3 cm x 24 cm <sup>(4)</sup> 04 cm x 20 cm <sup>(5)</sup>
Field of application	vaginal/cervical stump prolapse, concomitant cystocele/rectocele	
Surgical access	laparoscopic / open	
Surgical technique	colposacropexy / cervicosacropexy unilateral	
Fixation on vagina / cervix	sutures	
Fixation on sacrum	sutures / tacks	
Specially Warp-knitted Selvedges	●	
Shape stability	●	
Defined elasticity	●	
Visible technology	● (1,2,6)	● (3,4,5)
Polymer (monofilament)	PVDF	
Biocompatibility	●	
Ageing resistance	●	
Dynamometry	●	
Tear propagation resistance	●	
Classification (Klinge's classification [8])	1a	

### Application Examples:




Colpo-/cervicosacropexy

- unilateral



Colpo-/cervicosacropexy

- unilateral
- anterior/posterior mesh plasty  
(for concomitant  
cystocele/rectocele)



Colpo-/cervicosacropexy


- unilateral

Distributed by:

**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: info@zisanamed.ch

● Applies to all product sizes  
● Does not apply




Female Pelvic Organ Prolapse  
Vaginal/Cervical Stump Prolapse  
or Uterine Prolapse



**DynaMesh®-PRP** implants have been specially developed for pelvic floor reconstruction, in laparoscopic or open surgical technique, and serve to support and stabilise fascial structures, connective tissue and ligaments. The implants are used in the surgical treatment of a prolapse of the vaginal/cervical stump or uterine prolapse in the pectopexy technique.

## DynaMesh®-PRP

<b>DynaMesh®-PRP soft</b>	<sup>1)</sup> 03 cm x 15 cm	PV540315F1	BX = 1 piece
	<sup>1)</sup> 03 cm x 15 cm	PV540315F3	BX = 3 pieces
<b>DynaMesh®-PRP visible</b>	<sup>2)</sup> 03 cm x 15 cm	PV780315F1	BX = 1 piece
<b>DynaMesh®-PRP visible</b>	<sup>3)</sup> 03 cm x 18 cm	PV780318F1	BX = 1 piece
		PV780318F3	BX = 3 pieces
<b>DynaMesh®-PRP visible</b>	<sup>4)</sup> 17 cm x 15 cm	PV781715F1	BX = 1 piece
		PV781715F3	BX = 3 pieces




VI042xx	DynaMesh®-PRP - Animation: Pectopexy <a href="https://de.dyna-mesh.com/Vi042xx">https://de.dyna-mesh.com/Vi042xx</a>	
VI061xx	DynaMesh®-PRP - Animation: Hysteropectopexy - Anterior Fixation <a href="https://de.dyna-mesh.com/Vi061xx">https://de.dyna-mesh.com/Vi061xx</a>	
VI053xx	DynaMesh®-PRP - Animation: Hysteropectopexy - Posterior Fixation <a href="https://de.dyna-mesh.com/Vi053xx">https://de.dyna-mesh.com/Vi053xx</a>	
VI054xx	DynaMesh®-PRP - Animation: Pectopexy with Anterior & Posterior Mesh Repair <a href="https://de.dyna-mesh.com/Vi054xx">https://de.dyna-mesh.com/Vi054xx</a>	
VI069xx	DynaMesh® MRI - Animation: MRI Reconstruction with DynaMesh®-PRP visible <a href="https://de.dyna-mesh.com/Vi069xx">https://de.dyna-mesh.com/Vi069xx</a>	



## Use and Properties

Product	DynaMesh®-PRP soft <sup>(1)</sup> / visible <sup>(2)</sup> 03 cm x 15 cm	DynaMesh®-PRP visible 03 cm x 18 cm <sup>(3)</sup>	DynaMesh®-PRP visible 17 cm x 15 cm <sup>(4)</sup>
Field of application	vaginal/cervical stump or uterine prolapse	vaginal stump or uterine prolapse	vaginal/cervical stump prolapse, concomitant cystocele/rectocele
Surgical access	laparoscopic / open		
Surgical technique	pectopexy bilateral		
Fixation on vagina / cervix	sutures		
Fixation on pectineal ligament	sutures		
Specially Warp-knitted Selvedges	●		
Shape stability	●		
Defined elasticity	●		
Visible technology	● (1) / ● (2)		● (3,4)
Polymer (monofilament)	PVDF		
Biocompatibility	●		
Ageing resistance	●		
Dynamometry	●		
Tear propagation resistance	●		
Classification (Klinge's classification [8])	1a		




## Pectopexy Bilateral Fixation on the Pectineal Ligament

Fig. left:  
Apical mesh repair following hysterectomy with  
**DynaMesh®-PRP soft / visible (03 cm x 15 cm)**

Distributed by:


**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: info@zisanamed.ch

● Applies to all product sizes  
● Does not apply




Pectopexy after vaginal/cervical prolapse:

- Two implant sizes are available in the following dimensions  
**DynaMesh®-PRP soft / visible 03 cm x 15 cm** and  
**DynaMesh®-PRP visible 03 cm x 18 cm.**
- With greatly shortened vaginas, e.g., following a radical hysterectomy,  
**DynaMesh®-PRP visible 03 cm x 18 cm** can be optionally used.




Pectopexy after vaginal/cervical prolapse with concomitant cystocele and/or rectocele: (pulsion cystocele / rectocele)

- Additional stabilisation of the affected vaginal wall can be achieved with **DynaMesh®-PRP visible 17 cm x 15 cm.**



Pectopexy after uterine prolapse with uterine preservation:

- With a normal sized uterus, **DynaMesh®-PRP visible 03 cm x 18 cm** should be used and fixed in place on the posterior cervix.




Pectopexy after uterine prolapse with uterine preservation:

- With smaller uteri (below 100 g), anterior fixation of **DynaMesh®-PRP soft / visible 03 cm x 15 cm** can be selected as an alternative.

Distributed by:

**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: info@zisanamed.ch

Female Pelvic Organ Prolapse  
Cervical Stump Prolapse




DynaMesh®-CESA implants have been specially developed for pelvic floor reconstruction, and particularly for reinforcing or replacing the uterosacral ligaments, in laparoscopic or open surgical technique.

The implants are used in the treatment of a prolapse of the internal genitalia, such as a cervical stump prolapse.

## DynaMesh®-CESA


DynaMesh®-CESA	03 cm x 04 cm	PV740404F1	BX = 1 piece
		PV740404F3	BX = 3 pieces



### Use and Properties

Product	DynaMesh®-CESA
Field of application	cervical stump prolapse
Surgical access	laparoscopic / open
Surgical technique	cervicosacropexy (CESA) bilateral
Fixation on cervical stump	sutures
Fixation on sacrum	sutures / tacks
Specially Warp-knitted Selvedges	●
Shape stability	●
Defined elasticity	●
Visible technology	●
Polymer (monofilament)	PVDF
Biocompatibility	●
Ageing resistance	●
Dynamometry	●
Tear propagation resistance	●
Classification (Klinge's classification [8])	1a

● Applies to all product sizes




**DynaMesh®-CESA**  
(Cervico-SAcropexy)

The surgical technique CESA is a modified abdominal cervicosacropepy procedure (laparoscopic/open), in which the uterosacral ligaments are bilaterally reinforced or replaced by the implant.



**DynaMesh®-IVT02** instrument for **DynaMesh®-CESA** in retroperitoneal tape position through laparotomic access.  
Reusable instrument made of surgical steel.  
Length: 32 cm




- Extraperitoneal tunnelling
- Anatomically adapted to the pelvis
- Eyelet on instrument tip with slanted, atraumatic edges
- Use in laparoscopy
- Reusable instrument

VI094xx	DynaMesh®-CESA - Animation: Cervicosacropepy - Bilateral Fixation - Level Promontory <a href="https://de.dyna-mesh.com/Vi094xx">https://de.dyna-mesh.com/Vi094xx</a>	
VI084xx	DynaMesh®-CESA - Animation: Cervicosacropepy - Bilateral Fixation - Level S2 <a href="https://de.dyna-mesh.com/Vi084xx">https://de.dyna-mesh.com/Vi084xx</a>	

Distributed by:

**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: info@zisanamed.ch

Female Pelvic Organ Prolapse  
Vaginal Stump Prolapse




DynaMesh®-VASA implants have been specially developed for pelvic floor reconstruction, and particularly for reinforcing or replacing the uterosacral ligaments, in laparoscopic or open surgical technique.

The implants are used in the treatment of a prolapse of the internal genitalia, such as a vaginal stump prolapse.

## DynaMesh®-VASA


DynaMesh®-VASA	02 cm x 03 cm	PV740203F1	BX = 1 piece
		PV740203F3	BX = 3 pieces



### Use and Properties

Product	DynaMesh®-VASA
Field of application	vaginal stump prolapse
Surgical access	laparoscopic / open
Surgical technique	colposacropexy (VASA) bilateral
Fixation on vaginal stump	sutures
Fixation on sacrum	sutures / tacks
Specially Warp-knitted Selvedges	●
Shape stability	●
Defined elasticity	●
Visible technology	●
Polymer (monofilament)	PVDF
Biocompatibility	●
Ageing resistance	●
Dynamometry	●
Tear propagation resistance	●
Classification (Klinge's classification [8])	1a

● Applies to all product sizes




**DynaMesh®-VASA**  
(VAgino-SACropexy)

The surgical technique VASA is a modified abdominal colposacropexy procedure (laparoscopic/open), in which the uterosacral ligaments are bilaterally reinforced or replaced by the implant.



**DynaMesh®-IVT02** instrument for **DynaMesh®-VASA** in retroperitoneal tape position through laparotomic access.  
Reusable instrument made of surgical steel.  
Length: 32 cm




- Extraperitoneal tunnelling
- Anatomically adapted to the pelvis
- Eyelet on instrument tip with slanted, atraumatic edges
- Use in laparoscopy
- Reusable instrument

Distributed by:

**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: info@zisanamed.ch


Male Urinary Incontinence  
Stress Urinary Incontinence



DynaMesh®-PRM implants are used to support and stabilise connective tissue structures and ligaments. Common applications are suburethral slings for treating male stress urinary incontinence.

## DynaMesh®-PRM

DynaMesh®-PRM	04 cm x 03 cm	PV330453F1	BX = 1 piece
DynaMesh®-PRM visible	04 cm x 03 cm	PV730453F1	BX = 1 piece




### Use and Properties

Product	DynaMesh®-PRM	DynaMesh®-PRM visible
Field of application	stress urinary incontinence (SUI)	
Surgical access	perineal	
Surgical technique	Male Sling TOT - transobturator - outside-in	
Fixation	synthetic adhesives / sutures	
Specially Warp-knitted Selvedges		●
Shape stability		●
Defined elasticity		●
Visible technology	●	●
Polymer (monofilament)		PVDF
Biocompatibility		●
Ageing resistance		●
Dynamometry		●
Tear propagation resistance		●
Classification (Klinge's classification [8])		1a

● Applies to all product sizes  
● Does not apply

Application of the implant through perineal access  
Transobturator position



**DynaMesh®-IST03**  
Diameter: 5 cm

**DynaMesh®-IST02**  
Diameter: 7 cm

**DynaMesh®-IST03/-IST02:**

Instrument set consisting of two instruments (right and left side) for transobturator positioning using the outside-in technique.

Distributed by:

**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: info@zisanamed.ch



**Reusable Instruments**  
Manufactured from surgical steel (resterilisable)

For **transobturator** application


DynaMesh®-IST03

Surgical instrument

Diameter: 5 cm

IST03F1

BX = 1 set (l+r)




DynaMesh®-IST01

Surgical instrument

Diameter: 6 cm

IST01F1

BX = 1 set (l+r)




DynaMesh®-IST02

Surgical instrument

Diameter: 7 cm

IST02F1

BX = 1 set (l+r)




DynaMesh®-IVT01

Surgical instrument

IVT01F1

BX = 1 piece




For **retropubic** application

DynaMesh®-ISR01

Surgical instrument

ISR01F1

BX = 1 piece




For **laparotomical** application  
of DynaMesh®-CESA/-VASA

DynaMesh®-IVT02

Surgical instrument

IVT02F1

BX = 1 piece



Distributed by:

**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: info@zisanamed.ch

1. Klinge U, Klosterhalfen B, Ottinger AP, et al (2002)  
**PVDF as a new polymer for the construction of surgical meshes.**  
*Biomaterials* 23:3487-3493
2. Klink CD, Junge K, Binnebösel M, et al (2011)  
**Comparison of long-term biocompatibility of PVDF and PP meshes.**  
*J Invest Surg* 24:292-299. <https://doi.org/10.3109/08941939.2011.589883>
3. Gerullis H, Georgas E, Eimer C, et al (2011)  
**Evaluation of Biocompatibility of Alloplastic Materials: Development of a Tissue Culture In Vitro Test System.**  
*Surgical technology international* 21:21
4. Gerullis H, Klosterhalfen B, Borós M, et al (2013)  
**IDEAL in Meshes for Prolapse, Urinary Incontinence, and Hernia Repair.**  
*Surg Innov.* <https://doi.org/10.1177/1553350612472987>
5. Laroche G, Marois Y, Schwarz E, et al (1995)  
**Polyvinylidene fluoride monofilament sutures: can they be used safely for long-term anastomoses in the thoracic aorta?**  
*Artif Organs* 19:1190-1199
10. Berger D, Bientzle M (2008)  
**Polyvinylidene fluoride: a suitable mesh material for laparoscopic incisional and parastomal hernia repair!**  
*Hernia* 13:167-172. <https://doi.org/10.1007/s10029-008-0435-4>
11. Junge K, Binnebösel M, Rosch R, et al (2008)  
**Adhesion formation of a polyvinylidenfluoride/polypropylene mesh for intra-abdominal placement in a rodent animal model**  
*Surgical Endoscopy* 23:327-333. <https://doi.org/10.1007/s00464-008-9923-y>
20. Göretzlehner U, Müllen A (2007)  
**PVDF als Implantat-Werkstoff in der Urogynäkologie.**  
*BIOmaterialien* 8 (S1):28-29
27. Mary C, Marois Y, King MW, et al (1998)  
**Comparison of the in vivo behavior of polyvinylidene fluoride and polypropylene sutures used in vascular surgery.**  
*ASAIO J* 44:199-206
50. Roman S, Urbánková I, Callewaert G, et al (2016)  
**Evaluating Alternative Materials for the Treatment of Stress Urinary Incontinence and Pelvic Organ Prolapse: A Comparison of the In Vivo Response to Meshes Implanted in Rabbits.**  
*The Journal of Urology* 196:261-269. <https://doi.org/10.1016/j.juro.2016.02.067>
52. Silva RA, Silva PA, Carvalho ME (2007)  
**Degradation studies of some polymeric biomaterials: Polypropylene (PP) and polyvinylidene difluoride (PVDF).**  
*THERMEC 2006, Pts 1-5* 539-543:573-576
68. Conze J, Junge K, Weiss C, et al (2008)  
**New polymer for intra-abdominal meshes-PVDF copolymer.**  
*J Biomed Mater Res Part B Appl Biomater* 87:321-328. <https://doi.org/10.1002/jbm.b.31106>

## Implant Material

91. Hara T (2004)  
**Ten-Year Results of Anterior Chamber Fixation of the Posterior Chamber Intraocular Lens.**  
*Arch Ophthalmol* 122:1112. <https://doi.org/10.1001/archophth.122.8.1112>
93. Wang H, Klosterhalfen B, Müllen A, et al (2021)  
**Degradation resistance of PVDF mesh in vivo in comparison to PP mesh.**  
*J Mech Behav Biomed Mater* 119:104490. <https://doi.org/10.1016/j.jmbbm.2021.104490>
100. Karabulut A, Simavlı SA, Abban GM, et al (2016)  
**Tissue reaction to urogynecologic meshes: effect of steroid soaking in two different mesh models.**  
*Int Urogynecol J* 27:1583–1589. <https://doi.org/10.1007/s00192-016-3013-9>

6. Mühl T, Binnebösel M, Klinge U, Goedderz T (2008)  
**New objective measurement to characterize the porosity of textile implants.**  
*Journal of Biomedical Materials Research Part B: Applied Biomaterials* 84B:176–183. <https://doi.org/10.1002/jbm.b.30859>
8. Klinge U, Klosterhalfen B (2012)  
**Modified classification of surgical meshes for hernia repair based on the analyses of 1,000 explanted meshes.**  
*Hernia* 16:251–258. <https://doi.org/10.1007/s10029-012-0913-6>
25. Klosterhalfen B, Junge K, Klinge U (2005)  
**The lightweight and large porous mesh concept for hernia repair.**  
*Expert Rev Med Devices* 2:103–117. <https://doi.org/10.1586/17434440.2.1.103>
26. Otto J, Kaldenhoff E, Kirschner-Hermanns R, et al (2013)  
**Elongation of textile pelvic floor implants under load is related to complete loss of effective porosity, thereby favouring incorporation in scar plates.**  
*Journal of Biomedical Materials Research Part A n/a-n/a*. <https://doi.org/10.1002/jbm.a.34767>
38. Kaldenhoff E, Klinge U, Klosterhalfen B, et al (2013)  
**Von der Prolaps- zur Problempatientin: Schenken wir der Qualität von Netzimplantaten genügend Aufmerksamkeit?**  
*Der Gynäkologe* 46:469–476. <https://doi.org/10.1007/s00129-012-3124-4>
53. Zhu L-M, Schuster P, Klinge U (2015)  
**An overview of crucial mesh parameters.**  
*World Journal of Gastrointestinal Surgery*
102. Klinge U, Park J-K, Klosterhalfen B (2013)  
**The Ideal Mesh.**  
*Pathobiology* 80:169–175. <https://doi.org/10.1159/000348446>
103. Klosterhalfen B, Klinge U (2013)  
**Retrieval study at 623 human mesh explants made of polypropylene - impact of mesh class and indication for mesh removal on tissue reaction.**  
*Journal of Biomedical Materials Research Part B: Applied Biomaterials n/a-n/a*. <https://doi.org/10.1002/jbmb.32958>
104. Klinge U, Junge K, Spellerberg B, et al (2002)  
**Do multifilament alloplastic meshes increase the infection rate? Analysis of the polymeric surface, the bacteria adherence, and the in vivo consequences in a rat model.**  
*J Biomed Mater Res* 63:765–771. <https://doi.org/10.1002/jbm.10449>
105. Klinge U, Klosterhalfen B, Birkenhauer V, et al (2002)  
**Impact of polymer pore size on the interface scar formation in a rat model.**  
*J Surg Res* 103:208–214. <https://doi.org/10.1006/jsre.2002.6358>

## DynaMesh® visible

7. Hansen NL, Barabasch A, Distelmaier M, et al (2013)  
**First In-Human Magnetic Resonance Visualization of Surgical Mesh Implants for Inguinal Hernia Treatment.**  
*Invest Radiol.* <https://doi.org/10.1097/RLI.0b013e31829806ce>
29. Kuehnert N, Kraemer NA, Otto J, et al (2011)  
**In vivo MRI visualization of mesh shrinkage using surgical implants loaded with superparamagnetic iron oxides.**  
*Surgical Endoscopy* 26:1468-1475. <https://doi.org/10.1007/s00464-011-2057-7>
51. Köhler G, Pallwein-Prettner L, Lechner M, et al (2015)  
**First human magnetic resonance visualisation of prosthetics for laparoscopic large hiatal hernia repair.**  
*Hernia* 19:975-982. <https://doi.org/10.1007/s10029-015-1398-x>
54. Muysoms F, Beckers R, Kyle-Leinhase I (2018)  
**Prospective cohort study on mesh shrinkage measured with MRI after laparoscopic ventral hernia repair with an intraperitoneal iron oxide-loaded PVDF mesh.**  
*Surgical Endoscopy* 32:2822-2830. <https://doi.org/10.1007/s00464-017-5987-x>
56. Köhler G, Pallwein-Prettner L, Koch OO, et al (2015)  
**Magnetic Resonance-Visible Meshes for Laparoscopic Ventral Hernia Repair.**  
*JLS : Journal of the Society of Laparoendoscopic Surgeons* 19:e2014.00175. <https://doi.org/10.4293/JLS.2014.00175>
62. Köhler G, Wundsam H, Pallwein-Prettner L, et al (2015)  
**Magnetic resonance visible 3-D funnel meshes for laparoscopic parastomal hernia prevention and treatment.**  
*European Surgery* 47:127-132. <https://doi.org/10.1007/s10353-015-0319-7>
69. Kuehnert N, Otto J, Conze J, et al (2014)  
**Time-Dependent Changes of Magnetic Resonance Imaging-Visible Mesh Implants in Patients**
70. Hansen NL, Ciritsis A, Otto J, et al (2015)  
**Utility of Magnetic Resonance Imaging to Monitor Surgical Meshes: Correlating Imaging and Clinical Outcome of Patients Undergoing Inguinal Hernia Repair.**  
*Invest Radiol.* <https://doi.org/10.1097/RLI.0000000000000148>
71. Weyhe D, Klinge U, Uslar VN, et al (2019)  
**Follow Up Data of MRI-Visible Synthetic Meshes for Reinforcement in Large Hiatal Hernia in Comparison to None-Mesh Repair - A Prospective Cohort Study.**  
*Front Surg* 6:. <https://doi.org/10.3389/fsurg.2019.00017>
76. Lechner M, Meissnitzer M, Borhanian K, et al (2019)  
**Surgical and radiological behavior of MRI-depictable mesh implants after TAPP repair: the IRONMAN study.**  
*Hernia.* <https://doi.org/10.1007/s10029-019-02019-2>
90. Özveri E, Şanlı DET, Yıldırım D, et al (2020)  
**Magnetic resonance visualization of iron-loaded meshes in patients with pain after inguinal hernia repair.**  
*Hernia.* <https://doi.org/10.1007/s10029-020-02168-9>

## Implant Fixation

79. Villalobos RN, Mias MC, Gas C, et al (2019)  
**Atraumatic laparoscopic intraperitoneal mesh fixation using a new laparoscopic device: an animal experimental study.**  
*Hernia.* <https://doi.org/10.1007/s10029-019-02008-5>
86. Wilson P (2020)  
**Laparoscopic intraperitoneal onlay mesh (IPOM) repair using n-butyl-2-cyanoacrylate (Liquiband Fix8™) for mesh fixation: learning experience and short-medium term results.**  
*Hernia.* <https://doi.org/10.1007/s10029-020-02144-3>
97. Carus T (2021)  
**Die laparoskopische IPOM-Operation bei Nabel- und Bauchwandhernien – Netzfixierung in Klebetechnik.**  
6

## Female Pelvic Organ Prolapse

### DynaMesh®-PR / -PRR

57. Jan H, Ghai V, Thakar R (2018)  
**Simplified Laparoscopic Sacrohysteropexy.**  
*Journal of Minimally Invasive Gynecology* 25:1134. <https://doi.org/10.1016/j.jmig.2018.01.014>

## Female Pelvic Organ Prolapse

### DynaMesh®-PRS

42. Balsamo R, Illiano E, Zucchi A, et al (2018)

**Sacrocolpopexy with polyvinylidene fluoride mesh for pelvic organ prolapse:  
Mid term comparative outcomes with polypropylene mesh.**

*European Journal of Obstetrics & Gynecology and Reproductive Biology* 220:74–78. <https://doi.org/10.1016/j.ejogrb.2017.11.018>



## Female Pelvic Organ Prolapse

## DynaMesh®-PRP

22. Noé KG, Spüntrup C, Anapolski M (2013)  
**Laparoscopic pectopexy: a randomised comparative clinical trial of standard laparoscopic sacral colpo-cervicopexy to the new laparoscopic pectopexy. Short-term postoperative results.**  
*Arch Gynecol Obstet* 287:275–280. <https://doi.org/10.1007/s00404-012-2536-7>
32. Noé K-G, Schiermeier S, Alkatout I, Anapolski M (2015)  
**Laparoscopic Pectopexy: A Prospective, Randomized, Comparative Clinical Trial of Standard Laparoscopic Sacral Colpocervicopexy with the New Laparoscopic Pectopexy—Postoperative Results and Intermediate-Term Follow-Up in a Pilot Study.**  
*Journal of Endourology* 29:210–215. <https://doi.org/10.1089/end.2014.0413>
45. Kale A, Biler A, Terzi H, et al (2017)  
**Laparoscopic pectopexy: initial experience of single center with a new technique for apical prolapse surgery.**  
*International braz j urol* 43:903–909. <https://doi.org/10.1590/s1677-5538.ibju.2017.0070>
84. Noé GK, Schiermeier S, Papatthemelis T, et al (2020)  
**Prospective international multicenter pectopexy trial: Interim results and findings post surgery.**  
*European Journal of Obstetrics & Gynecology and Reproductive Biology* 244:81–86. <https://doi.org/10.1016/j.ejogrb.2019.10.022>
98. Noé GK, Schiermeier S, Papatthemelis T, et al (2021)  
**Prospective International Multicenter Pelvic Floor Study: Short-Term Follow-Up and Clinical Findings for Combined Pectopexy and Native Tissue Repair.**  
*JCM* 10:217. <https://doi.org/10.3390/jcm10020217>

## Female Pelvic Organ Prolapse

### DynaMesh®-CESA / -VASA

33. Jaeger et al. (2016)  
**Does the Patients Age have an Influence on the Outcome of Cesa (Cervico-Sacropexy) and Vasa (Vagino-Sacropexy) for the Treatment of Urinary Incontinence in Women?**  
*| Open Access | OMICS International*
34. Rajshekhar S, Mukhopadhyay S, Morris E (2016)  
**Early safety and efficacy outcomes of a novel technique of sacrocolpopexy for the treatment of apical prolapse.**  
*International Journal of Gynecology and Obstetrics 0: <https://doi.org/10.1016/j.ijgo.2016.05.007>*
35. Ludwig S, Stumm M, Mallmann P, Jager W (2016)  
**Surgical replacement of the uterosacral and pubourethral ligaments as treatment for urgency urinary incontinence.**  
*Austin J Womens Health 3:1019*
36. Joukhadar R, Meyberg-Solomayer G, Hamza A, et al (2015)  
**A Novel Operative Procedure for Pelvic Organ Prolapse Utilizing a MRI-Visible Mesh Implant: Safety and Outcome of Modified Laparoscopic Bilateral Sacropexy.**  
*BioMed Research International 2015:1-9. <https://doi.org/10.1155/2015/860784>*
73. Cassis C, Mukhopadhyay S, Morris E (2019)  
**Standardizing abdominal sacrocolpopexy for the treatment of apical prolapse: One year on.**  
*Int J Gynecol Obstet ijgo.12935. <https://doi.org/10.1002/ijgo.12935>*
85. Rexhepi S, Rexhepi E, Stumm M, et al (2018)  
**Laparoscopic Bilateral Cervicosacropexy and Vaginosacropexy: New Surgical Treatment Option in Women with Pelvic Organ Prolapse and Urinary Incontinence.**  
*Journal of Endourology 32:1058-1064. <https://doi.org/10.1089/end.2018.0474>*

## Female Urinary Incontinence

21. Klinge U, Binneboesel M, Kuschel S, Schuessler B (2007)  
**Demands and properties of alloplastic implants for the treatment of stress urinary incontinence.**  
*Expert Rev Med Devices* 4:349–359. <https://doi.org/10.1586/17434440.4.3.349>
31. Naumann G, Albrich S, Skala C, et al (2012)  
**Single-Incision Slings (SIS) - a New Option for the Surgical Treatment of Female Stress Urinary Incontinence.**  
*Geburtshilfe und Frauenheilkunde* 72:125–131. <https://doi.org/10.1055/s-0031-1298275>
39. Ludwig S, Stumm M, Mallmann P, Jäger W (2016)  
**TOT 8/4: A Way to Standardize the Surgical Procedure of a Transobturator Tape.**  
*BioMed Research International* 2016:1–4. <https://doi.org/10.1155/2016/4941304>
40. Najjari L, Hennemann J, Kirschner-Hermanns R, et al (2014)  
**Visualization of Polypropylene and Polyvinylidene Fluoride Slings in Perineal Ultrasound and Correlation with Clinical Outcome.**  
*BioMed Research International*. <https://doi.org/10.1155/2014/181035>
41. Sabadell J, Larrain F, Gracia-Perez-Bonfils A, et al (2016)  
**Comparative study of polyvinylidene fluoride and polypropylene suburethral slings in the treatment of female stress urinary incontinence: PVDF/polypropylene in suburethral slings.**  
*Journal of Obstetrics and Gynaecology Research* 42:291–296. <https://doi.org/10.1111/jog.12899>
44. Ludwig S, Stumm M (2016)  
**Surgical Treatment of Urgency Urinary Incontinence, OAB (Wet), Mixed Urinary Incontinence, and Total Incontinence by Cervicosacropexy or Vaginosacropexy.**  
*Gynecology & Obstetrics* 6:. <https://doi.org/10.4172/2161-0932.1000404>
49. Najjari L, Gräf CM, Kupec T, et al (2016)  
**Tomographic Ultrasound Imaging to Control the Placement of Tension-Free Transobturator Tape in Female Urinary Stress Incontinence.**  
*BioMed Research International* 2016:1–6. <https://doi.org/10.1155/2016/6495858>
83. Ludwig S, Becker I, Mallmann P, Jäger W (2019)  
**Comparison of Solifenacin and Bilateral Apical Fixation in the Treatment of Mixed and Urgency Urinary Incontinence in Women: URGE 1 Study, A Randomized Clinical Trial.**  
*In Vivo* 33:1949–1957. <https://doi.org/10.21873/invivo.11690>
99. Sabadell J, Pereda-Núñez A, Ojeda-de-los-Santos F, et al (2021)  
**Polypropylene and polyvinylidene fluoride transobturator slings for the treatment of female stress urinary incontinence: 1-Year outcomes from a multicentre randomized trial.**  
*Neurourology and Urodynamics* 40:475–482. <https://doi.org/10.1002/nau.24586>

## Male Urinary Incontinence

81. Costa Cruz DSL da, D´Ancona CAL, Silva Filho WP da, et al (2020)  
**Parameters of 2-Dimensional Perineal Ultrasonography Before and After Male Sling Procedure for Urinary Incontinence After Radical Prostatectomy.**  
*Urology 136:257–262. <https://doi.org/10.1016/j.urology.2019.10.004>*

## Internal test-reports

TR1. CF\_F02 **“Biocompatibility”**

TR10. I1\_F03-01\_SIS1\_X **Stability and Elongation** (*bench test*)

TR11. I1\_F03-02\_SIS1\_X **Porosity and Formstability** (*bench test*)

TR12. I1\_F03-03\_SIS1\_X **Formstability** (*bench test*)

**Some of our products may not be available in your country. Please contact your local distributor for more information.**






**[www.dyna-mesh.com](http://www.dyna-mesh.com)**

hergestellt durch / manufactured by /  
fabriqué par / fabricado por / fabbricato da  
FEG Textiltechnik  
Forschungs- und Entwicklungsgesellschaft mbH  
Prager Ring 70  
52070 Aachen, Germany



Distributed by:

**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: [info@zisanamed.ch](mailto:info@zisanamed.ch)



**Download  
Catalogue**




**DynaMesh**<sup>®</sup>

Implants for the surgical treatment of  
Female Urinary Incontinence  
Female Pelvic Organ Prolapse  
Male Urinary Incontinence



made  
in  
Germany

Tailored Implants  
Made of **PVDF**



hergestellt durch / manufactured by: FEG Textiltechnik Forschungs- und Entwicklungsgesellschaft mbH  
Prager Ring 70, 52070 Aachen, Germany · [www.dyna-mesh.com](http://www.dyna-mesh.com)

Distributed by:

**Zi+Sanamed GmbH**  
Bernstrasse 18  
CH-2555 Brügg  
phone +41 793310767  
email: [info@zisanamed.ch](mailto:info@zisanamed.ch)


**LiSTO.academy** is a comprehensive platform providing excellent customised surgical education and trainings. Through cooperations with experienced surgeons and renowned experts from around the world, **LiSTO.academy** enables physicians to achieve best patient outcomes using DynaMesh® implants.



**IN-PERSON**



**VIRTUAL**



**HYBRID**

Visit us online:  
[www.listo.academy](http://www.listo.academy)

