
# DynaMesh<sup>®</sup>

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## Implantat Material

1. Klinge U, Klosterhalfen B, Ottinger AP, et al (2002)  
**PVDF as a new polymer for the construction of surgical meshes.**  
*Biomaterials* 23:3487-3493
2. Klink CD, Junge K, Binnebösel M, et al (2011)  
**Comparison of long-term biocompatibility of PVDF and PP meshes.**  
*J Invest Surg* 24:292-299. <https://doi.org/10.3109/08941939.2011.589883>
3. Gerullis H, Georgas E, Eimer C, et al (2011)  
**Evaluation of Biocompatibility of Alloplastic Materials: Development of a Tissue Culture In Vitro Test System.**  
*Surgical technology international* 21:21
4. Gerullis H, Klosterhalfen B, Borós M, et al (2013)  
**IDEAL in Meshes for Prolapse, Urinary Incontinence, and Hernia Repair.**  
*Surg Innov.* <https://doi.org/10.1177/1553350612472987>
5. Laroche G, Marois Y, Schwarz E, et al (1995)  
**Polyvinylidene fluoride monofilament sutures: can they be used safely for long-term anastomoses in the thoracic aorta?**  
*Artif Organs* 19:1190-1199
10. Berger D, Bientzle M (2008)  
**Polyvinylidene fluoride: a suitable mesh material for laparoscopic incisional and parastomal hernia repair!**  
*Hernia* 13:167-172. <https://doi.org/10.1007/s10029-008-0435-4>
11. Junge K, Binnebösel M, Rosch R, et al (2008)  
**Adhesion formation of a polyvinylidenfluoride/polypropylene mesh for intra-abdominal placement in a rodent animal model**  
*Surgical Endoscopy* 23:327-333. <https://doi.org/10.1007/s00464-008-9923-y>
20. Göretzlehner U, Müllen A (2007)  
**PVDF als Implantat-Werkstoff in der Urogynäkologie.**  
*BIOmaterialien* 8 (S1):28-29
27. Mary C, Marois Y, King MW, et al (1998)  
**Comparison of the in vivo behavior of polyvinylidene fluoride and polypropylene sutures used in vascular surgery.**  
*ASAIO J* 44:199-206
50. Roman S, Urbánková I, Callewaert G, et al (2016)  
**Evaluating Alternative Materials for the Treatment of Stress Urinary Incontinence and Pelvic Organ Prolapse: A Comparison of the In Vivo Response to Meshes Implanted in Rabbits.**  
*The Journal of Urology* 196:261-269. <https://doi.org/10.1016/j.juro.2016.02.067>
52. Silva RA, Silva PA, Carvalho ME (2007)  
**Degradation studies of some polymeric biomaterials: Polypropylene (PP) and polyvinylidene difluoride (PVDF).**  
*THERMEC 2006, Pts 1-5* 539-543:573-576
68. Conze J, Junge K, Weiss C, et al (2008)  
**New polymer for intra-abdominal meshes-PVDF copolymer.**  
*J Biomed Mater Res Part B Appl Biomater* 87:321-328. <https://doi.org/10.1002/jbm.b.31106>

## Implantat Material

91. Hara T (2004)  
**Ten-Year Results of Anterior Chamber Fixation of the Posterior Chamber Intraocular Lens.**  
*Arch Ophthalmol* 122:1112. <https://doi.org/10.1001/archophth.122.8.1112>
93. Wang H, Klosterhalfen B, Müllen A, et al (2021)  
**Degradation resistance of PVDF mesh in vivo in comparison to PP mesh.**  
*J Mech Behav Biomed Mater* 119:104490. <https://doi.org/10.1016/j.jmbbm.2021.104490>
100. Karabulut A, Simavlı SA, Abban GM, et al (2016)  
**Tissue reaction to urogynecologic meshes: effect of steroid soaking in two different mesh models.**  
*Int Urogynecol J* 27:1583–1589. <https://doi.org/10.1007/s00192-016-3013-9>

6. Mühl T, Binnebösel M, Klinge U, Goedderz T (2008)  
**New objective measurement to characterize the porosity of textile implants.**  
*Journal of Biomedical Materials Research Part B: Applied Biomaterials* 84B:176–183. <https://doi.org/10.1002/jbm.b.30859>
8. Klinge U, Klosterhalfen B (2012)  
**Modified classification of surgical meshes for hernia repair based on the analyses of 1,000 explanted meshes.**  
*Hernia* 16:251–258. <https://doi.org/10.1007/s10029-012-0913-6>
25. Klosterhalfen B, Junge K, Klinge U (2005)  
**The lightweight and large porous mesh concept for hernia repair.**  
*Expert Rev Med Devices* 2:103–117. <https://doi.org/10.1586/17434440.2.1.103>
26. Otto J, Kaldenhoff E, Kirschner-Hermanns R, et al (2013)  
**Elongation of textile pelvic floor implants under load is related to complete loss of effective porosity, thereby favouring incorporation in scar plates.**  
*Journal of Biomedical Materials Research Part A n/a-n/a*. <https://doi.org/10.1002/jbm.a.34767>
38. Kaldenhoff E, Klinge U, Klosterhalfen B, et al (2013)  
**Von der Prolaps- zur Problempatientin: Schenken wir der Qualität von Netzimplantaten genügend Aufmerksamkeit?**  
*Der Gynäkologe* 46:469–476. <https://doi.org/10.1007/s00129-012-3124-4>
53. Zhu L-M, Schuster P, Klinge U (2015)  
**An overview of crucial mesh parameters.**  
*World Journal of Gastrointestinal Surgery*
102. Klinge U, Park J-K, Klosterhalfen B (2013)  
**The Ideal Mesh.**  
*Pathobiology* 80:169–175. <https://doi.org/10.1159/000348446>
103. Klosterhalfen B, Klinge U (2013)  
**Retrieval study at 623 human mesh explants made of polypropylene - impact of mesh class and indication for mesh removal on tissue reaction.**  
*Journal of Biomedical Materials Research Part B: Applied Biomaterials n/a-n/a*. <https://doi.org/10.1002/jbmb.32958>
104. Klinge U, Junge K, Spellerberg B, et al (2002)  
**Do multifilament alloplastic meshes increase the infection rate? Analysis of the polymeric surface, the bacteria adherence, and the in vivo consequences in a rat model.**  
*J Biomed Mater Res* 63:765–771. <https://doi.org/10.1002/jbm.10449>
105. Klinge U, Klosterhalfen B, Birkenhauer V, et al (2002)  
**Impact of polymer pore size on the interface scar formation in a rat model.**  
*J Surg Res* 103:208–214. <https://doi.org/10.1006/jsre.2002.6358>

## DynaMesh® visible

7. Hansen NL, Barabasch A, Distelmaier M, et al (2013)  
**First In-Human Magnetic Resonance Visualization of Surgical Mesh Implants for Inguinal Hernia Treatment.**  
*Invest Radiol.* <https://doi.org/10.1097/RLI.0b013e31829806ce>
29. Kuehnert N, Kraemer NA, Otto J, et al (2011)  
**In vivo MRI visualization of mesh shrinkage using surgical implants loaded with superparamagnetic iron oxides.**  
*Surgical Endoscopy* 26:1468-1475. <https://doi.org/10.1007/s00464-011-2057-7>
51. Köhler G, Pallwein-Prettner L, Lechner M, et al (2015)  
**First human magnetic resonance visualisation of prosthetics for laparoscopic large hiatal hernia repair.**  
*Hernia* 19:975-982. <https://doi.org/10.1007/s10029-015-1398-x>
54. Muysoms F, Beckers R, Kyle-Leinhase I (2018)  
**Prospective cohort study on mesh shrinkage measured with MRI after laparoscopic ventral hernia repair with an intraperitoneal iron oxide-loaded PVDF mesh.**  
*Surgical Endoscopy* 32:2822-2830. <https://doi.org/10.1007/s00464-017-5987-x>
56. Köhler G, Pallwein-Prettner L, Koch OO, et al (2015)  
**Magnetic Resonance-Visible Meshes for Laparoscopic Ventral Hernia Repair.**  
*JLS : Journal of the Society of Laparoendoscopic Surgeons* 19:e2014.00175. <https://doi.org/10.4293/JLS.2014.00175>
62. Köhler G, Wundsam H, Pallwein-Prettner L, et al (2015)  
**Magnetic resonance visible 3-D funnel meshes for laparoscopic parastomal hernia prevention and treatment.**  
*European Surgery* 47:127-132. <https://doi.org/10.1007/s10353-015-0319-7>
69. Kuehnert N, Otto J, Conze J, et al (2014)  
**Time-Dependent Changes of Magnetic Resonance Imaging-Visible Mesh Implants in Patients**
70. Hansen NL, Ciritsis A, Otto J, et al (2015)  
**Utility of Magnetic Resonance Imaging to Monitor Surgical Meshes: Correlating Imaging and Clinical Outcome of Patients Undergoing Inguinal Hernia Repair.**  
*Invest Radiol.* <https://doi.org/10.1097/RLI.000000000000148>
71. Weyhe D, Klinge U, Uslar VN, et al (2019)  
**Follow Up Data of MRI-Visible Synthetic Meshes for Reinforcement in Large Hiatal Hernia in Comparison to None-Mesh Repair - A Prospective Cohort Study.**  
*Front Surg* 6:. <https://doi.org/10.3389/fsurg.2019.00017>
76. Lechner M, Meissnitzer M, Borhanian K, et al (2019)  
**Surgical and radiological behavior of MRI-depictable mesh implants after TAPP repair: the IRONMAN study.**  
*Hernia.* <https://doi.org/10.1007/s10029-019-02019-2>
90. Özveri E, Şanlı DET, Yıldırım D, et al (2020)  
**Magnetic resonance visualization of iron-loaded meshes in patients with pain after inguinal hernia repair.**  
*Hernia.* <https://doi.org/10.1007/s10029-020-02168-9>

## Implantat Fixation

79. Villalobos RN, Mias MC, Gas C, et al (2019)  
**Atraumatic laparoscopic intraperitoneal mesh fixation using a new laparoscopic device: an animal experimental study.**  
*Hernia.* <https://doi.org/10.1007/s10029-019-02008-5>
86. Wilson P (2020)  
**Laparoscopic intraperitoneal onlay mesh (IPOM) repair using n-butyl-2-cyanoacrylate (Liquiband Fix8™) for mesh fixation: learning experience and short-medium term results.**  
*Hernia.* <https://doi.org/10.1007/s10029-020-02144-3>
97. Carus T (2021)  
**Die laparoskopische IPOM-Operation bei Nabel- und Bauchwandhernien – Netzfixierung in Klebetechnik.**  
6

## Leistenhernie

16. Junge K, Binnebösel M, Kauffmann C, et al (2010)  
**Damage to the spermatic cord by the Lichtenstein and TAPP procedures in a pig model.**  
*Surgical Endoscopy* 25:146-152. <https://doi.org/10.1007/s00464-010-1148-1>
67. Garcia-Pastor P, Porrero-Carro J, et al. (2018)  
**Prospective Multicenter Blinded Randomized Study Comparing PP and PVDF Mesh Implants in Lichtenstein Procedure with Respect to Pain and Recurrence.**  
*JSM Surgical Procedures* 1:
72. Guadalajara Jurado JF, Suárez Grau JM, Bellido Luque JA, et al (2016)  
**Initial experience in laparoscopic bilateral inguinal hernia repair (TEP) with new anatomical mesh with large pore and low weight (Dynamesh Endolap) in short stay (6 months follow-up).**  
*Ambulatory Surgery* 22:
90. Özveri E, Şanlı DET, Yıldırım D, et al (2020)  
**Magnetic resonance visualization of iron-loaded meshes in patients with pain after inguinal hernia repair.**  
*Hernia*. <https://doi.org/10.1007/s10029-020-02168-9>
94. Ramser M, Baur J, Keller N, et al (2021)  
**Robotische Hernienchirurgie: Teil I: Robotische Leistenhernienversorgung (r TAPP). Videobeitrag und Ergebnisse einer Kohortenstudie an 302 operierten Hernien.**  
*Chirurg*. <https://doi.org/10.1007/s00104-021-01425-6>
101. The HerniaSurge Group (2018)  
**International guidelines for groin hernia management.**  
*Hernia*. <https://doi.org/10.1007/s10029-017-1668-x>

## Reparatur der Bauchwandhernie

10. Berger D, Bientzle M (2008)  
**Polyvinylidene fluoride: a suitable mesh material for laparoscopic incisional and parastomal hernia repair!**  
*Hernia* 13:167-172. <https://doi.org/10.1007/s10029-008-0435-4>
14. Berger D, Bientzle M (2006)  
**Principles of laparoscopic repair of ventral hernias.**  
*European Surgery* 38:393-398. <https://doi.org/10.1007/s10353-006-0284-2>
54. Muysoms F, Beckers R, Kyle-Leinhase I (2018)  
**Prospective cohort study on mesh shrinkage measured with MRI after laparoscopic ventral hernia repair with an intraperitoneal iron oxide-loaded PVDF mesh.**  
*Surgical Endoscopy* 32:2822-2830. <https://doi.org/10.1007/s00464-017-5987-x>
56. Köhler G, Pallwein-Prettner L, Koch OO, et al (2015)  
**Magnetic Resonance-Visible Meshes for Laparoscopic Ventral Hernia Repair.**  
*JLS : Journal of the Society of Laparoendoscopic Surgeons* 19:e2014.00175. <https://doi.org/10.4293/JLS.2014.00175>
58. Verbo A, Pafundi P, Manno A, et al (2016)  
**Polyvinylidene Fluoride Mesh (PVDF, DynaMesh®-IPOM) in The Laparoscopic Treatment of Incisional Hernia: A Prospective Comparative Trial versus Gore® ePTFE DUALMESH® Plus.**  
*Surgical technology international* 28:147-151
87. Sánchez-Arteaga A, Tallón-Aguilar L, Tinoco-González J, et al (2020)  
**Use of polyvinylidene fluoride (PVDF) meshes for ventral hernia repair in emergency surgery.**  
*Hernia*. <https://doi.org/10.1007/s10029-020-02209-3>

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55. Kohler A, Lavanchy JL, Lenoir U, et al (2019)  
**Effectiveness of Prophylactic Intraperitoneal Mesh Implantation for Prevention of Incisional Hernia in Patients Undergoing Open Abdominal Surgery: A Randomized Clinical Trial.**  
*JAMA Surgery* 154:109. <https://doi.org/10.1001/jamasurg.2018.4221>
74. Bravo-Salva A, González-Castillo AM, Vela-Polanco FF, et al (2019)  
**Incidence of Incisional Hernia After Emergency Subcostal Unilateral Laparotomy: Does Augmentation Prophylaxis Play a Role?**  
*World J Surg.* <https://doi.org/10.1007/s00268-019-05282-7>
89. Pereira JA, Pera M, López-Cano M, et al (2019)  
**Hernias at the Extraction Incision After Laparoscopic Colon and Rectal Resection: Influence of Incision Location and Use of Prophylactic Mesh.**  
*Cirugía Española (English Edition)* 97:20–26

## Reparatur der Parastomalhernie

9. Berger D, Bientzle M (2007)  
**Laparoscopic Repair of Parastomal Hernias: A Single Surgeon's Experience in 66 Patients.**  
*Diseases of the Colon & Rectum* 50:1668–1673. <https://doi.org/10.1007/s10350-007-9028-z>
12. Berger D (2010)  
**Laparoskopische Reparatur der parastomalen Hernie.**  
*Der Chirurg* 81:988–992. <https://doi.org/10.1007/s00104-010-1933-3>
60. Fischer I, Wundsam H, Mitteregger M, Köhler G (2017)  
**Parastomal Hernia Repair with a 3D Funnel Intraperitoneal Mesh Device and Same-Sided Stoma Relocation: Results of 56 Cases.**  
*World Journal of Surgery* 41:3212–3217. <https://doi.org/10.1007/s00268-017-4130-4>
62. Köhler G, Wundsam H, Pallwein-Prettner L, et al (2015)  
**Magnetic resonance visible 3-D funnel meshes for laparoscopic parastomal hernia prevention and treatment.**  
*European Surgery* 47:127–132. <https://doi.org/10.1007/s10353-015-0319-7>
63. Köhler G, Emmanuel K (2017)  
**Laparoscopic stoma relocation for parastomal hernia treatment by using a magnetic resonance visible three-dimensional implant.**  
*ANZ Journal of Surgery* 87:411–412. <https://doi.org/10.1111/ans.12899>
64. Köhler G, Fischer I, Wundsam H (2018)  
**A Novel Technique for Parastomal Hernia Repair Combining a Laparoscopic and Ostomy-Opening Approach.**  
*Journal of Laparoendoscopic & Advanced Surgical Techniques* 28:209–214.  
<https://doi.org/10.1089/lap.2017.0313>
65. Köhler G, Mayer F, Wundsam H, et al (2015)  
**Changes in the Surgical Management of Parastomal Hernias Over 15 Years: Results of 135 Cases.**  
*World Journal of Surgery* 39:2795–2804. <https://doi.org/10.1007/s00268-015-3187-1>
66. Zhang H, Xie J, Miao J, Wu H (2016)  
**Hybrid Approaches for Complex Parastomal Hernia Repair.**  
*Journal of the College of Physicians and Surgeons Pakistan* 26:72–73
75. Köhler G (2019)  
**Prinzipien und Parallelen der Prävention und Reparatur parastomaler Hernien mit Netzen.**  
*Chirurg*. <https://doi.org/10.1007/s00104-019-01047-z>
77. Szczepkowski M, Skoneczny P, Przywózka A, et al (2015)  
**Leading article: Methods paper New minimally invasive technique of parastomal hernia repair - methods and review.**  
*wiitm* 1:1–7. <https://doi.org/10.5114/wiitm.2015.50052>
78. Tully KH, Roghmann F, Pastor J, et al (2019)  
**Parastomal Hernia Repair With 3-D Mesh Implants After Radical Cystectomy and Ileal Conduit Urinary Diversion - A Single-center Experience Using a Purpose Made Alloplastic Mesh Implant.**  
*Urology* 131:245–249. <https://doi.org/10.1016/j.urology.2019.05.006>

## Hernien

## Reparatur der Parastomalhernie

92. Bustos-Jiménez M, Martín-Cartes JA (2020)  
**Surgical Treatment of Parastomal Hernias by Using A 3D Mesh.**  
*05:6*
96. Cartes JAM, Bustos-Jiménez M, Tamayo-López MJ (2018)  
**Parastomal Hernia: A More and More Frequent Surgical Challenge.**  
*General Surgery 3:5*

## Prävention der Parastomalhernie

15. Berger D (2007)  
**Prevention of parastomal hernias by prophylactic use of a specially designed intraperitoneal onlay mesh (Dynamesh IPST®).**  
*Hernia* 12:243–246. <https://doi.org/10.1007/s10029-007-0318-0>
59. Conde-Muñoz R, Díez J-L, Martínez A, et al (2017)  
**Preventing parastomal hernias with systematic intraperitoneal specifically designed mesh.**  
*BMC Surgery* 17:. <https://doi.org/10.1186/s12893-017-0237-7>
61. Köhler G, Hofmann A, Lechner M, et al (2016)  
**Prevention of parastomal hernias with 3D funnel meshes in intraperitoneal onlay position by placement during initial stoma formation.**  
*Hernia* 20:151–159. <https://doi.org/10.1007/s10029-015-1380-7>
62. Köhler G, Wundsam H, Pallwein-Prettner L, et al (2015)  
**Magnetic resonance visible 3-D funnel meshes for laparoscopic parastomal hernia prevention and treatment.**  
*European Surgery* 47:127–132. <https://doi.org/10.1007/s10353-015-0319-7>
75. Köhler G (2019)  
**Prinzipien und Parallelen der Prävention und Reparatur parastomaler Hernien mit Netzen.**  
*Chirurg.* <https://doi.org/10.1007/s00104-019-01047-z>
80. López-Borao J, Madrazo-González Z, Kreisler E, Biondo S (2019)  
**Prevention of parastomal hernia after abdominoperineal excision with a prophylactic three-dimensional funnel mesh.**  
*Colorectal Dis* 21:1326–1334. <https://doi.org/10.1111/codi.14738>
88. Mäkäräinen-Uhlbäck EJ, Klintrup KHB, Vierimaa MT, et al (2020)  
**Prospective, Randomized Study on the Use of Prosthetic Mesh to Prevent a Parastomal Hernia in a Permanent Colostomy: Results of a Long-term Follow-up.**  
*Diseases of the Colon & Rectum* 63:678–684. <https://doi.org/10.1097/DCR.0000000000001599>
95. Ammann Y, Widmann B, Sparr M, et al (2021)  
**Prophylactic Funnel Mesh to Prevent Parastomal Hernia in Permanent End Colostomy: A Retrospective Cohort Study.**  
*Colorectal Dis.* <https://doi.org/10.1111/codi.15817>

## Reparatur der Hiatushernie

51. Köhler G, Pallwein-Prettner L, Lechner M, et al (2015)  
**First human magnetic resonance visualisation of prosthetics for laparoscopic large hiatal hernia repair.**  
*Hernia* 19:975-982. <https://doi.org/10.1007/s10029-015-1398-x>
71. Weyhe D, Klinge U, Uslar VN, et al (2019)  
**Follow Up Data of MRI-Visible Synthetic Meshes for Reinforcement in Large Hiatal Hernia in Comparison to None-Mesh Repair – A Prospective Cohort Study.**  
*Front Surg* 6:. <https://doi.org/10.3389/fsurg.2019.00017>

## Weibliche Beckenbodensenkung

### DynaMesh®-PR / -PRR

57. Jan H, Ghai V, Thakar R (2018)

**Simplified Laparoscopic Sacrohysteropexy.**

*Journal of Minimally Invasive Gynecology* 25:1134. <https://doi.org/10.1016/j.jmig.2018.01.014>

## Weibliche Beckenbodensenkung

### DynaMesh®-PRS

42. Balsamo R, Illiano E, Zucchi A, et al (2018)

**Sacrocolpopexy with polyvinylidene fluoride mesh for pelvic organ prolapse:  
Mid term comparative outcomes with polypropylene mesh.**

*European Journal of Obstetrics & Gynecology and Reproductive Biology* 220:74–78. <https://doi.org/10.1016/j.ejogrb.2017.11.018>

## Weibliche Beckenbodensenkung

## DynaMesh®-PRP

22. Noé KG, Spüntrup C, Anapolski M (2013)  
**Laparoscopic pectopexy: a randomised comparative clinical trial of standard laparoscopic sacral colpo-cervicopexy to the new laparoscopic pectopexy. Short-term postoperative results.**  
*Arch Gynecol Obstet* 287:275–280. <https://doi.org/10.1007/s00404-012-2536-7>
32. Noé K-G, Schiermeier S, Alkatout I, Anapolski M (2015)  
**Laparoscopic Pectopexy: A Prospective, Randomized, Comparative Clinical Trial of Standard Laparoscopic Sacral Colpocervicopexy with the New Laparoscopic Pectopexy—Postoperative Results and Intermediate-Term Follow-Up in a Pilot Study.**  
*Journal of Endourology* 29:210–215. <https://doi.org/10.1089/end.2014.0413>
45. Kale A, Biler A, Terzi H, et al (2017)  
**Laparoscopic pectopexy: initial experience of single center with a new technique for apical prolapse surgery.**  
*International braz j urol* 43:903–909. <https://doi.org/10.1590/s1677-5538.ibju.2017.0070>
84. Noé GK, Schiermeier S, Papatthemelis T, et al (2020)  
**Prospective international multicenter pectopexy trial: Interim results and findings post surgery.**  
*European Journal of Obstetrics & Gynecology and Reproductive Biology* 244:81–86. <https://doi.org/10.1016/j.ejogrb.2019.10.022>
98. Noé GK, Schiermeier S, Papatthemelis T, et al (2021)  
**Prospective International Multicenter Pelvic Floor Study: Short-Term Follow-Up and Clinical Findings for Combined Pectopexy and Native Tissue Repair.**  
*JCM* 10:217. <https://doi.org/10.3390/jcm10020217>

## Weibliche Beckenbodensenkung

## DynaMesh®-CESA / -VASA

33. Jaeger et al. (2016)  
**Does the Patients Age have an Influence on the Outcome of Cesa (Cervico-Sacropexy) and Vasa (Vagino-Sacropexy) for the Treatment of Urinary Incontinence in Women?**  
*| Open Access | OMICS International*
34. Rajshekhar S, Mukhopadhyay S, Morris E (2016)  
**Early safety and efficacy outcomes of a novel technique of sacrocolpopexy for the treatment of apical prolapse.**  
*International Journal of Gynecology and Obstetrics 0: <https://doi.org/10.1016/j.ijgo.2016.05.007>*
35. Ludwig S, Stumm M, Mallmann P, Jager W (2016)  
**Surgical replacement of the uterosacral and pubourethral ligaments as treatment for urgency urinary incontinence.**  
*Austin J Womens Health 3:1019*
36. Joukhadar R, Meyberg-Solomayer G, Hamza A, et al (2015)  
**A Novel Operative Procedure for Pelvic Organ Prolapse Utilizing a MRI-Visible Mesh Implant: Safety and Outcome of Modified Laparoscopic Bilateral Sacropexy.**  
*BioMed Research International 2015:1-9. <https://doi.org/10.1155/2015/860784>*
73. Cassis C, Mukhopadhyay S, Morris E (2019)  
**Standardizing abdominal sacrocolpopexy for the treatment of apical prolapse: One year on.**  
*Int J Gynecol Obstet ijgo.12935. <https://doi.org/10.1002/ijgo.12935>*
85. Rexhepi S, Rexhepi E, Stumm M, et al (2018)  
**Laparoscopic Bilateral Cervicosacropexy and Vaginosacropexy: New Surgical Treatment Option in Women with Pelvic Organ Prolapse and Urinary Incontinence.**  
*Journal of Endourology 32:1058-1064. <https://doi.org/10.1089/end.2018.0474>*

## Weibliche Harninkontinenz

21. Klinge U, Binneboesel M, Kuschel S, Schuessler B (2007)  
**Demands and properties of alloplastic implants for the treatment of stress urinary incontinence.**  
*Expert Rev Med Devices* 4:349–359. <https://doi.org/10.1586/17434440.4.3.349>
31. Naumann G, Albrich S, Skala C, et al (2012)  
**Single-Incision Slings (SIS) - a New Option for the Surgical Treatment of Female Stress Urinary Incontinence.**  
*Geburtshilfe und Frauenheilkunde* 72:125–131. <https://doi.org/10.1055/s-0031-1298275>
39. Ludwig S, Stumm M, Mallmann P, Jäger W (2016)  
**TOT 8/4: A Way to Standardize the Surgical Procedure of a Transobturator Tape.**  
*BioMed Research International* 2016:1–4. <https://doi.org/10.1155/2016/4941304>
40. Najjari L, Hennemann J, Kirschner-Hermanns R, et al (2014)  
**Visualization of Polypropylene and Polyvinylidene Fluoride Slings in Perineal Ultrasound and Correlation with Clinical Outcome.**  
*BioMed Research International*. <https://doi.org/10.1155/2014/181035>
41. Sabadell J, Larrain F, Gracia-Perez-Bonfils A, et al (2016)  
**Comparative study of polyvinylidene fluoride and polypropylene suburethral slings in the treatment of female stress urinary incontinence: PVDF/polypropylene in suburethral slings.**  
*Journal of Obstetrics and Gynaecology Research* 42:291–296. <https://doi.org/10.1111/jog.12899>
44. Ludwig S, Stumm M (2016)  
**Surgical Treatment of Urgency Urinary Incontinence, OAB (Wet), Mixed Urinary Incontinence, and Total Incontinence by Cervicosacropexy or Vaginosacropexy.**  
*Gynecology & Obstetrics* 6:. <https://doi.org/10.4172/2161-0932.1000404>
49. Najjari L, Gräf CM, Kupec T, et al (2016)  
**Tomographic Ultrasound Imaging to Control the Placement of Tension-Free Transobturator Tape in Female Urinary Stress Incontinence.**  
*BioMed Research International* 2016:1–6. <https://doi.org/10.1155/2016/6495858>
83. Ludwig S, Becker I, Mallmann P, Jäger W (2019)  
**Comparison of Solifenacin and Bilateral Apical Fixation in the Treatment of Mixed and Urgency Urinary Incontinence in Women: URGE 1 Study, A Randomized Clinical Trial.**  
*In Vivo* 33:1949–1957. <https://doi.org/10.21873/invivo.11690>
99. Sabadell J, Pereda-Núñez A, Ojeda-de-los-Santos F, et al (2021)  
**Polypropylene and polyvinylidene fluoride transobturator slings for the treatment of female stress urinary incontinence: 1-Year outcomes from a multicentre randomized trial.**  
*Neurourology and Urodynamics* 40:475–482. <https://doi.org/10.1002/nau.24586>

## Männliche Harninkontinenz

81. Costa Cruz DSL da, D´Ancona CAL, Silva Filho WP da, et al (2020)  
**Parameters of 2-Dimensional Perineal Ultrasonography Before and After Male Sling Procedure for Urinary Incontinence After Radical Prostatectomy.**  
*Urology 136:257–262. <https://doi.org/10.1016/j.urology.2019.10.004>*

## Internal test-reports

- TR1. CF\_F02 **“Biocompatibility”**
- TR10. I1\_F03-01\_SIS1\_X **Stability and Elongation** (*bench test*)
- TR11. I1\_F03-02\_SIS1\_X **Porosity and Formstability** (*bench test*)
- TR12. I1\_F03-03\_SIS1\_X **Formstability** (*bench test*)

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Dahlhausen GmbH